
 DIAGNOSTIC SERVICES OF MANITOBA SERVICES DE DIAGNOSTIC DU MANITOBA	Equipment Malfunction and Corrective Action Record: Storage Equipment and Plasma Thawer		Document # F160-QCFORM-15
	Approved by: 	Effective Date: 31-MAR-2011	Version # 01 Source Document: Manitoba Transfusion Quality Manual for Blood Banks

Equipment _____

Model # _____ Serial # _____ Location: _____

Date: _____ Time: _____ Reported by: _____

Part A: Equipment Malfunction Identification and Actions

If this is a satellite refrigerator, did the nursing unit call you? Yes No

State Problem	Temperatures		
	Chart	Digital Controller (Internal Thermometer)	Independent Thermometer
<input type="checkbox"/> Refrigerator alarm triggered. Temperature at time of alarm			
<input type="checkbox"/> Freezer alarm triggered. Temperature at time of alarm			
<input type="checkbox"/> Platelet Incubator alarm triggered. Temperature at time of alarm			
<input type="checkbox"/> Room alarm triggered. Temperature at time of alarm			
<input type="checkbox"/> Plasma thawer outside of 30-37°C range. Temperature at time of alarm			
<input type="checkbox"/> Malfunction of Continuous Recording Graph only <ul style="list-style-type: none"> <input type="checkbox"/> Document temperatures every four hours until resolved 			
<input type="checkbox"/> Other			

If applicable: Length of time out of acceptable range _____ Remote Alarm Activated: Yes No

Cause of alarm: Door ajar Other (Specify): _____

Describe action(s) taken. Include temperatures taken at 10 or 15 minute intervals:

Action	Time	Tech	Temperatures		
			Chart	Digital Controller Internal Thermometer	Independent Thermometer

Was problem resolved? Yes No

If No:
Document on refrigerator as "STOP Service required" as per policy 100-10-07 Appendix 4
Proceed to Part B.



Equipment Malfunction and Corrective Action Record: Storage Equipment and Plasma Thawer

Document # F160-QCFORM-15

Version # 01

Effective Date: 31-MAR-2011

Part B: Malfunction (specify) _____

Interim Corrective Actions

Relocation to on site alternate storage: No **If Yes = Complete Chart below:**

Location Details	Alternate Storage Equipment Name	Model#	Serial #	Location
Date (dd/mmm/yy)	Time Product Removed	Relocation Completion Time	Product Type(s)	Signature

Alternate storage temperature monitoring implemented (check all that apply)

<input type="checkbox"/> 4 hour check required (No continuous chart recorder available – use alternative Storage Temperature Record Form)	<input type="checkbox"/> Continuous monitoring available (4 hour not required)
<input type="checkbox"/> Daily temperature Record	<input type="checkbox"/> Capacity for alarm activation and monitoring (required to meet back-up criteria.

If No:

Inter-facility transfer activated (policy *Inter-facility Shipping of Blood, Blood Components and Derivatives*) Yes

Site: _____

Inter-facility documentation completed Yes

Attach all supporting interim corrective action documents.

Part C: Equipment Return to Service

Validation Performed as per policy *Validation of Blood, Blood Components and Derivatives temperature monitored storage equipment* and applicable form for documentation Yes

Placement of Product in equipment Date: _____ Initials: _____

Attach all supporting return to service documents.

Part D: Supervisory Review

Was the service company or maintenance called? Yes No

Was the problem resolved? Yes No

Corrective action taken: _____

Maintenance documents attached to report? Yes No

Reviewed by (signature/date): _____